


How to Install Satellite Licenses


Please click [here](#) for system requirements.

For help, contact our support team in [Europe](#), [Americas](#), and [APAC](#).

- 1) Access the Customer Portal [here](#). Select an appropriate Site if you have more than one.
- 2) Download the latest version of Imaris (9.5).
 - a. Select 'Download'.
 - b. Choose the appropriate file format.



- 3) The following options are only available to the Site Leader. Set the duration of the Satellite License.
 - a. Select 'Manager User'.
 - b. Set the duration (2 to 30 days. Default option = 7 days).
 - c. Save.



4) Generate Satellite Licenses.

- a. Click on 'Satellite Licensing'.
- b. Click on the link at the bottom of the page to create a new satellite license.

Won Yung Choi
US-BPI.Internal

HOME
APPLICATIONS
PRODUCTS
CONTACT
EVENTS
TRAINING

PORTAL
User Settings
Satellite Licensing
Report a Bug
Manage User
Product Status
Maintenance
News Letters
Transfer Rights
Download
Issue Voucher
Training Archive
Memory
logout

Satellite Licensing

Please read changes to Satellite License due SARS-CoV-2 [here](#).

This service is provided for those individuals who need to show the use of the Imaplane software as a part of a presentation or trade show while traveling or need to be able to use the software at home or away from the office for a small amount of time.

We reserve the right to monitor the use of these temporary licenses.

The Imaplane team is happy to provide alternative temporary licenses for microscopy and image analysis courses. Please contact your local [sales and support teams](#) to discuss how we can support your course and to make licensing arrangements.

Once you have created the license, you will receive an identification code called OrderID and is formatted like "abcd-1234-abcd-1234". The OrderID will activate Imaplane on the new computer. A guide on how to install Imaplane can be found [here](#).

To create a new satellite license, please click [here](#).

HostID	Article	OrderID	Issued	Start Date	Expiration
--------	---------	---------	--------	------------	------------

- c. Select the module or package, and fill out the form (see an example below).
- d. Click on "Create Satellite Licenses"
- e. **Note the OrderID below and copy it to save it in a text file.** It will also be listed on the Satellite Licensing at the bottom of the screen above (after step b) if you accidentally closed the window.
- f. **MAKE SURE TO click on the last 'here' button.**

Won Yung Choi
US-BPI.Internal

HOME
APPLICATIONS
PRODUCTS
CONTACT
EVENTS
TRAINING

PORTAL
User Settings
Satellite Licensing
Report a Bug
Manage User
Product Status
Maintenance
News Letters
Transfer Rights
Download
Issue Voucher
Training Archive
Memory
logout

BITPLANE AG
SCIENTIFIC SOLUTIONS

Satellite Licensing

Please select product you like to create a satellite license:

Articles with Maintenance Contract

☐ ClearView Deconvolution 9.3 (0 of 2 licenses used.)
☐ ClearView Deconvolution 9.5 (1 of 2 licenses used.)
☐ Imaris 9.5 Single Full (Old with AA) (2 of 2 licenses used. If you need more satellite license, please contact [sales and support](#).)
☒ Imaris Stitcher 9.5 (1 of 2 licenses used.)

Articles without Maintenance Contract

☐ AutoAligner 6.0 (0 of 1 licenses used.)
☐ ClearView Deconvolution 9.5 (1 of 2 licenses used.)
☐ Imaris 8.4 Single Full (SF) (0 of 1 licenses used.)
☐ Imaris 9.2 Single Full (Old with AA) (0 of 1 licenses used.)
☐ Imaris 9.5 Single Full (SF) (0 of 3 licenses used.)
☐ Imaris Stitcher 9.2 (0 of 1 licenses used.)
☐ Imaris Stitcher 9.5 (1 of 3 licenses used.)

Please state your reason for generating a satellite license: *

To access Imaris from the home office.

Which statement describes the best your purpose for creating a license?

d

Won Yung Choi
US-BPI.Internal

HOME
APPLICATIONS
PRODUCTS
CONTACT
EVENTS
TRAINING

PORTAL
User Settings
Satellite Licensing
Report a Bug
Manage User
Product Status
Maintenance
News Letters
Transfer Rights
Download
Issue Voucher
Training Archive
Memory
logout

BITPLANE AG
SCIENTIFIC SOLUTIONS


Satellite Licensing

Satellite license was created successfully!


The OrderID, which registers this satellite license, is: **y4gx-78dh-j23x**
 At any time you can come back to the portal and look it up here again.

f Click [here](#) to return to the overview!

- 4) Launch Imaris. You need to have Administrator access.
a. Click on Configurator.



b. Click on 'Add License Location'



c. Select 'Download Node-Locked License'.

d. Paste Order ID acquired earlier in the field below under "Step 2". Click "Check Order-ID" to the right – you may have to be patient as this step may take a bit.

e. Select all the relevant modules under "New Licenses".

f. Click "OK". Then "Finish". Start Imaris and Enjoy your Satellite License!

